



Non-precious PFM

We are a Chinese dental lab. We provide custom-made dental service for dental institutions at home and abroad, and provide oem service for other dental factories. PFM is one of our main products. Made from materials certified by FDA and the European Union. PFM consists of non-precious metal porcelain, semi-precious metal porcelain and precious metal porcelain. Non-precious PFM.

China Non-precious PFM Outsourcing Dental Lab, Manufacturer, Non-precious PFM Lab to Lab Service



Description

Porcelain fused to metal crowns & bridges (PFM crowns) can be referred to as full-cast crowns & bridges which has porcelain fused on most or all parts of the tooth. ... Metal alloys with a high melting temperature are used to fuse porcelain to the surface and preventing the metal from melting.

Porcelain fused to metal crowns & bridges are more affordable than all ceramic. Even though they are not as aesthetic as porcelain crowns they are still a good candidate for back teeth restoration.



Introduction: PFM-Porcelain Fused to Metal Nickel Free

Material: Argeloy N.P. Supreme Metal(Product code:116422) /Vita vmk95

Shenzhen Wanmei Dental Laboratory Co.,LTD.

Tel: +86-755-29760054

E-mail: info@wmdentallab.com

Related products



[PFM Crowns & Bridge](#)



[Full Cast Crown & Bridge](#)



[All Ceramic](#)



[eMax](#)



[Custom Abutments](#)



[Uper Restoration](#)



[Flexible Denture](#)



[Metal Framework](#)

Material List

Material List

Type	Name of material	CE No.	Lot No.	Manufacturer	Composition
Non-precious Alloy	Argeloy N.P.Supreme	CE0086	14034132	Argen(USA)	Co60.%,Cr27%,W5%,Mo6%
Precious Alloy	Argedent 74	CE0086	718846622	Argen(USA)	Au74%,In2.5%,Pt4.0%,Pd7.6%
	Argedent 52	CE0086	71977557	Argen(USA)	Au52.5%,Pd26.9%,Ag16%,In2.5%
Semi-precious	Argistar 43+	CE0086	7545643	Argen(USA)	Au2%,Pd42.9%,Ag46%,In8%
Metal Framework	Argeloy N.P.Partial	CE0086	13869812	Argen(USA)	Co64%,Cr28%,Mo6%,
	Vitalium 2000	CE01663	N001830	Dentsply(USA)	Co63.8%,Cr28%,Mo5.9%
Acrylic Resin	TCS Unbreakable	CE0470	73021314	TCS(USA)	Acrylic Resin
	Vertex	CE120	XT016907	Castapress	Acrylic Resin
	Lucitone 199	CE PR017600	131114	Dentsply(USA)	Acrylic Resin
All Ceramic	E max Press 3	CE0123	P19900	Ivoclar(Liechtenstein)	SiO279%,LiO15%,K2O5%,P2O5%,ZrO2%,ZrO3%
	Zirconia	HD601143840001	W1802260NG-1	Aidite(China)	ZrO2,Y2O3,AL2O3
	Full Contour Zirconia	HD601143840001	W171023QGC-1	Aidite(China)	ZrO2,Y2O3,AL2O3
Ceramic Powder	E max Ceram 3	CE0123	P62370	Ivoclar(Liechtenstein)	SiO265%,Al2O38%,Na2O6%,ZrO3%
	Vita Ceram	CE0124	LOT45190	Vita(Germany)	SiO265%,Al2O315%

The business keeps to the operation concept "scientific management, premium quality and efficiency primacy, customer supreme for Good quality Non-precious PFM From Chinese Dental, We have now experienced manufacturing facilities with a lot more than 100 workforce. So we could Shenzhen Wanmei Dental Laboratory Co.,LTD.

guarantee short lead time and top quality assurance.

Good quality China Dental Implants Product, Non-precious PFM, Selling our goods causes no risks and brings high returns to your company instead. It is our consistent pursuit to create value for clients. Our company is looking for agents sincerely. What are you waiting for? Come and join us. Now or never.